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# Samsung Health Monitor

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**The Product described in below is medical device**

# Samsung Health Monitor

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### (BP Monitor App)

129, Samsung-ro,  
Yeongtong-gu,  
Suwon-si, Gyeonggi-do,  
Republic of Korea

**Samsung Electronic Co., Ltd.**

[www.samsung.com](http://www.samsung.com)

Please read this instruction  
for use carefully to use the  
product properly

# Instruction for Use

2020. 05. 25 (Rev. A)

## 1 INTENDED USE

The BP Monitor App is intended for over-the-counter (OTC) use in adults 22 years and older. It may not provide accurate results in pregnant women. The BP data displayed by the BP Monitor App is intended for informational use only. The user is not intended to interpret or take clinical action based on the device output without consultation with a qualified healthcare professional and validation with a standard measurement method. The BP Monitor App is not intended to replace traditional methods for diagnosis or treatment of hypertension.

## 2 PRODUCT OVERVIEW

The Samsung Blood Pressure Monitor Application (BP Monitor App) is a software-only, mobile medical application intended for use with a compatible Samsung Galaxy Watch and Phone to create, record, store and display blood pressure information. The BP Monitor App determines the systolic and diastolic blood pressure as well as pulse rate and provides values and history charts following calibration with a blood pressure cuff. These data are only captured when the user is at rest. The BP Monitor App uses signals obtained from an optical sensor that measures signals non-invasively from the wrist.

## 3 WARNINGS

Please follow the warnings listed below.

|   |  |
|---|--|
|  | <p>The BP Monitor App cannot diagnose hypertension, other conditions, or check for signs of a heart attack. The BP Monitor App is not meant to replace traditional methods of diagnosis or treatment by a qualified healthcare professional. If you think you are having a medical emergency, contact your local emergency services immediately.</p> |
|  | <p>DO NOT change your medications or dosage without first consulting your doctor.</p>  |
|  | <p>Results displayed by the BP Monitor App are for <b>informational purposes only</b>. DO NOT interpret or take clinical action based on the BP Monitor App measurements without first consulting with a qualified healthcare professional.</p>  |
|  | <p>DO NOT use the BP Monitor App if you are pregnant.</p>  |
|  | <p>DO NOT use the BP Monitor App if you have any of the following:</p> <ul style="list-style-type: none"> <li>• Arrhythmia</li> <li>• Prior heart failure or heart attack</li> <li>• Peripheral vascular disease or compromised circulation</li> <li>• Valvular disease (diseases concerning the aortic valve)</li> <li>• Cardiomyopathy</li> </ul>  |

|   |   |
|---|---|
|   | <ul style="list-style-type: none"> <li>• Other known cardiovascular disease</li> <li>• End-stage Renal Disease (ESRD)</li> <li>• Diabetes</li> <li>• Neurological disorder (for example, tremors)</li> <li>• Clotting disorder or you are taking prescribed blood thinners</li> <li>• Dark tattoo(s) on the wrist where you will wear the Galaxy Watch</li> </ul> |
|  | Samsung BP Monitor app analyzes signal acquired using optic-based sensor (a photoplethysmography sensor, PPG). Therefore, the signal may be affected by light absorption by the skin, the amount of blood flowing in blood vessels, and contamination of the sensor.  |
|  | The Samsung Blood Pressure App calibrates and measures blood pressure based on the blood pressure value of the arm wearing the cuff. If difference of blood pressure values between both arms is large* (more than 10mmHg), please consult a health care professional for further evaluation  |

**CAUTIONS**

Please follow the precautions below. If you do not follow them, the BP Monitor App may not be able to measure blood pressure accurately.

|   |  |
|---|--|
|  | DO NOT use a Galaxy Watch for BP measurements if the Galaxy Watch has been calibrated for use by another person.   |
|  | DO NOT take BP measurements during any physical activity.  |
|  | DO NOT use the BP Monitor App if you are younger than 22 years old.  |
|  | DO NOT take BP measurements when the Galaxy Watch is close to strong electromagnetic fields (for example, magnetic resonance imaging (MRI) or X-Ray equipment, electromagnetic anti-theft systems, and metal detectors). |
|  | For security reasons, always pair the Galaxy phone and Galaxy watch via Bluetooth in a private (home-based) setting. It is NOT recommended that pairing be done in a public space.                                       |
|  | DO NOT take BP measurements during a medical procedure (for example, surgery or external defibrillation procedures).   |
|  | DO NOT take BP measurements when the Galaxy Watch is outside of the following operating conditions: <ul style="list-style-type: none"> <li>• Temperature: 54 °F - 104 °F / 12 °C - 40 °C</li> </ul>                      |

|   |  |
|---|--|
|   | <ul style="list-style-type: none"> <li>Humidity: 30% – 90% relative humidity</li> </ul>  |
|  | <p>The BP Monitor App's range for blood-pressure readings is as follows:</p> <ul style="list-style-type: none"> <li>Calibration Range - Systolic: 80–170 mmHg, Diastolic: 50–110 mmHg</li> <li>Measurement Range - Systolic: 70–180 mmHg, Diastolic: 40–120 mm Hg</li> </ul> |

## 4 Getting Started

Before you can use the BP Monitor App, you must complete pairing your phone with a Samsung Galaxy watch:

1. Set your watch in 'Get Started mode'.
2. Install and run Galaxy Wearable App on your Galaxy phone
3. Register your Galaxy watch and pair with each other

After pairing Galaxy phone and watch, you must do the following:

4. Set-up the BP Monitor App on your Galaxy Watch
5. Install the Samsung Health Monitor App on your Galaxy phone
6. Calibrate your Galaxy Watch with an MFDS(Ministry of Food and Drug Safety)-approved(certified) upper-arm, cuff-based blood pressure monitor

The BP Monitor App requires a compatible Galaxy Watch paired with a compatible Galaxy phone.

### 4.1 Setting up the BP Monitor App

Before you can start measuring your blood pressure using the BP Monitor App, you must set up your BP Monitor App profile and calibrate your Galaxy Watch.

1. Ensure that your Galaxy Watch is paired via Bluetooth with the Galaxy phone that you will use for your blood pressure measurements.
2. If the BP Monitor App is installed on your Galaxy Watch, open the BP Monitor App and follow the onscreen set-up instructions.

If the BP App is not already installed on your Galaxy phone, check that your Galaxy Watch is a Galaxy Watch Active2 model or one of more recent Samsung Galaxy Watch models. For more information about Galaxy Watch Active2, check the website, <https://www.samsung.com/sec/wearables/all-wearables/?galaxy-watch+galaxy-watch-active>. If your Galaxy Watch is compatible, upgrade to the latest Maintenance Release (MR) update using the Galaxy Wearable application on your Galaxy phone. After you have installed the update, open the Samsung Health Monitor App and follow the onscreen set-up instructions.

► In order to set up the BP Monitor App Profile on Your Galaxy Phone

1. Download the Samsung Health Monitor App from the Galaxy Store.
2. Complete and confirm the profile set-up and agree to the Terms of Use.
3. Select the Blood Pressure option and calibrate your Galaxy Watch.

## 4.2 Calibrating Your Galaxy Watch

To ensure more accurate measurements, be sure to calibrate your Galaxy Watch with an MFDS-approved (or certified) upper-arm, cuff-based blood pressure monitor. Before you start the calibration, read the instructions that accompany the cuff-based blood pressure monitor and make sure you know how to use it.

**Note 1:** You need the MFDS-approved (or certified) upper-arm, cuff-based blood pressure monitor only to calibrate your Galaxy Watch. After you calibrate your Galaxy Watch, you do not need to use the cuff-based blood pressure monitor when you use the BP Monitor App on your Galaxy Watch to take your blood pressure until recalibration.

**Note 2:** Calibration on the Galaxy Watch customizes the BP Monitor App specifically—and only—for the person wearing the paired Galaxy Watch during calibration.

To ensure more accurate calibration, do the following:

- 30 minutes before doing the calibration, do not take any alcohol or caffeine, and do not smoke, exercise, or bathe.
- Make sure your arm and wrist are dry and free of excessive perspiration or skin lotion.
- Do the calibration indoors in a quiet place.
- Put a comfortable chair next to a table.
- Sit in the chair with your back well-supported, your legs uncrossed, and your feet flat on the ground.
- Let your hands and forearms rest on the table.
- Rest in this position for at least 5 minutes before starting the calibration.
- When you start the calibration, stay still and do not move your arm or talk while the cuff-based blood pressure monitor performs the measurement.
- During the measurement, breathe normally. Do not try to deepen or slow down your breath.

To calibrate your Galaxy Watch:

1. Place the blood pressure cuff on an upper arm.
2. Wear the Samsung Galaxy Watch on the wrist of **the other** arm. Make sure the Galaxy Watch strap is snug around your wrist, but not too tight.
3. Place your Galaxy phone on the table, within easy reach.
4. On your Galaxy phone, open the BP Monitor App and follow the onscreen instructions to do the following:
  - ① Start the blood pressure measurement on the cuff-based blood pressure monitor.
  - ② Measurement on the Galaxy Watch will begin automatically.
  - ③ Enter the cuff-based blood pressure monitor reading in the phone BP Monitor App.
  - ④ Repeat Steps ① - ③ two more times (for a total of three measurements) to finish the Galaxy Watch calibration.

**Note:** Make sure that you finish the three measurements required for calibration continuously. To keep your Galaxy Watch properly calibrated, you must calibrate it every 28 days, as directed by the BP Monitor App. You may also calibrate your Galaxy Watch any time before then by selecting **Recalibrate** from the menu on the Blood Pressure History screen. If the wearer of the Galaxy Watch is changed or the wearing position is changed to the other arm, the watch must be recalibrated.

## 5 Measurement

### 5.1 Measuring Your Blood Pressure Using Your Galaxy Watch

After you have calibrated your Galaxy Watch with the MFDS- approved (or certified) upper-arm, cuff-based blood pressure monitor, you are ready to use the BP Monitor App on your Galaxy Watch to measure your blood pressure. You will not need to use the cuff-based monitor again until the next calibration in 28 days.

**Note:** After you calibrate your Galaxy Watch, you do not need to use the cuff-based monitor when you use the BP Monitor App on your Galaxy Watch to take your blood pressure until recalibration.

To ensure a more accurate BP measurement, do the following:

- Wear your Galaxy Watch on the same wrist as during the calibration and make sure the Galaxy Watch strap is snug around your wrist, but not too tight.
- 30 minutes before doing the measurement, do not have any alcohol or caffeine, and do not smoke, exercise, or bathe.
- Make sure your arm and wrist are free of excessive perspiration or skin lotion.
- Do the measurement indoors in a quiet place.
- Put a comfortable chair next to a table.
- Sit in the chair with your back well-supported, your legs uncrossed, and your feet flat on the ground.
- Let your hands and forearms rest on the table.
- Rest in this position for at least 5 minutes before starting the measurement.
- When you start the measurement, stay still and do not move your arm or talk while the Galaxy Watch takes the measurement.
- During the measurement, breathe normally. Do not try to deepen or slow down your breath.

▶ To measure your blood pressure

1. Press the **Measure** button on the Galaxy Watch.
2. If the measurement is successful, BP measurements including systolic and diastolic readings appear on the Galaxy Watch face. See “Reviewing Your Blood Pressure Measurements” below.  
If the measurement is unsuccessful (for example, because you moved or because the signal was not strong enough), the Galaxy Watch displays an error message.

## 6 Reviewing Your Blood Pressure Measurements

You can review your blood pressure measurements on your Galaxy Watch and on your Galaxy phone.

### 6.1 On Your Galaxy Watch

After the BP Monitor App on the Galaxy Watch successfully measures your blood pressure, the Galaxy Watch displays the following results:

- Systolic blood pressure (Systolic BP)
- Diastolic blood pressure (Diastolic BP)
- Pulse Rate (heart rate)

The table below offers a general guideline for blood pressure.

Blood Pressure Classification\* in adults (age 22 and above)

| BP Category                    | Systolic BP  |     | Diastolic BP |
|--------------------------------|--------------|-----|--------------|
| Normal                         | <120 mmHg    | and | <80 mmHg     |
| Elevated                       | 120-129 mmHg | and | <80 mmHg     |
| Prehypertension                | 130-139 mmHg | or  | 80-89 mmHg   |
| Hypertension Stage 1           | 140-159mmHg  | or  | 90-99mmHg    |
| Hypertension Stage 2           | >=160mmHg    | or  | >=100mmHg    |
| Isolated Systolic Hypertension | >=140mmHg    | and | <90mmHg      |

\* Individuals with Systolic BP and Diastolic BP in 2 different categories should be designated to the higher BP category. (e.g. If Sys.=110mmHg (Normal) and Dia.=88mmHg (Prehypertension), It should be designated as Prehypertension)

\* Source: Korean Society of Hypertension, 2018 High Blood Pressure Clinical Practice Guideline

## WARNING

Do not change your medications or dosage based on the numbers provided by the BP App. Always talk to your doctor first.

### 6.2 On Your Galaxy Phone

When a Bluetooth connection is established between your Galaxy Watch and the Galaxy phone you used to calibrate the Galaxy Watch, your measurement results are synced to the Samsung Health Monitor App on the Galaxy phone. You can share this information with your doctor for review or consultation.

► To view your blood pressure measurements on your phone

Open the Samsung Health Monitor App and swipe up the BP Monitor Card to see the history chart

**Note:** The number of results you can store on your Galaxy phone depends on the size of available storage memory in your Galaxy phone.

## 7 SAFETY AND PERFORMANCE

The Samsung BP Monitor App's ability to accurately measure blood pressure was tested in a clinical trial with 85 subjects. Manual blood pressure cuff measurements were taken by two medical professionals and compared to the simultaneously measured blood pressure values by the BP Monitor App. This was done following a modified protocol derived from the draft blood pressure testing standard ISO 81060-2:2018. The Samsung BP Monitor App measurements were shown to provide performance sufficient for the device's intended use. Specifically, the performance of the Phase 1 (Accuracy), Phase 2 (Induced Changes), and Phase 3 (BP Stability) tests yielded the results shown in the table below.

| Type                | Performance  | Test Protocol   |
|---------------------|--|---|
| Range               | Systolic BP : 70~180 mmHg<br>Diastolic BP : 40~120 mmHg<br>Pulse rate : 31 ~ 179bpm  |   |
| BP accuracy         | [Differences of Means]<br>Systolic BP : Average $\pm$ 5 mmHg under<br>Diastolic BP : Average $\pm$ 5 mmHg under<br>[Standard deviation]<br>Systolic BP : 8 mmHg under<br>Diastolic BP : 8 mmHg under | After calibrating with an auscultatory method, BP accuracy was assessed according to ISO 81060-2.   |
| Pulse rate accuracy | Under $\pm$ 5 %  | Pulse rate accuracy was assessed using a simulator, a similar method to the validation method of MFDS-approved (or certified) automatic blood pressure monitor. |

## 8 TROUBLESHOOTING

If you have problems using your BP Monitor App, check the list below for possible solutions.

**Problem:** I cannot find the BP Monitor App on my Galaxy Watch.

**Solution:**

- Check that your Galaxy Watch is a Galaxy Watch Active2 or more recent Galaxy Watch model that supports the BP Monitor App.
- If your Galaxy Watch is the right one, make sure it is upgraded to the latest Maintenance Release (MR) using the Galaxy Wearable application.

**Problem:** I don't see an option for measuring my blood pressure.

**Solution:**

- Make sure that you set up your profile on the paired-phone Samsung Health Monitor App. See "To set up the BP Monitor App Profile on Your Galaxy Phone" in this document.
- Make sure that you successfully have calibrated the Galaxy Watch with an MFDS-approved (or certified) upper-arm, cuff-based BP monitor. See "Calibrating your Galaxy Watch" on this document.
- If you have calibrated your Galaxy Watch, make sure that the calibration is up to date. The Galaxy Watch must be recalibrated every 28 days. Check the BP Monitor App screen of either Galaxy Watch or Galaxy phone to see the calibration status.

**Problem:** I cannot get a successful blood pressure reading.

**Solution:**

- Make sure the Galaxy Watch is worn snugly on your wrist (not loose or too tight). The bottom surface of the Galaxy Watch should touch the skin but should not be uncomfortable.
- Make sure you sit still and relaxed, with both feet flat on the ground, during the blood pressure measurement.
- Make sure you do not move your arm or talk during the blood pressure measurement.
- Rest the arm with the Galaxy Watch on a table while you take the measurement.
- Make sure your arm and wrist are free of excessive perspiration or skin lotion.
- If your hands and arms are cold, remove the Galaxy Watch and rub the wrist gently to warm up the skin. Put the Galaxy Watch back on and try again to take a measurement.
- During the measurement, breathe comfortably. Do not deliberately take deep breaths or slow breaths.
- If you can't make successful measurement repeatedly, please contact Samsung Service Center(phone: 1588-3366 or <https://www.samsungsvc.co.kr>)

**Problem:** I am getting an abnormally high or low BP reading.

**Solution:**

- If you have unusual symptoms or think you are having a medical emergency, contact your doctor or local emergency services immediately.
- If you do not think you are having a medical emergency, take at least two more blood pressure measurements and follow the guidance below before you take the measurements:
  - Do not have caffeine or alcohol 30 minutes before you start the measurement.
  - Do not smoke, exercise, or bathe 30 minutes before you start the measurement.
  - Rest quietly for at least 5 minutes before you start the measurement.
  - Stay rested during the measurement.
- If you are still getting abnormal readings, contact your doctor or local emergency services immediately for guidance.

**Problem:** I cannot sync my BP data from the Galaxy Watch to the Galaxy phone BP Monitor Card.

**Solution:**

- Make sure your Galaxy Watch is paired, using Bluetooth, with the Galaxy phone via the Galaxy Wearable application.
- Make sure you are using the same Galaxy phone that you used to set up your profile and calibrate the Galaxy Watch.
- Make sure you have at least 100 MBs of storage available on the Galaxy phone.